

# OUR DEVICES IN COMPARISON

Features	Teena	Daysy	Lady-Comp
Cycle overview	with phases	with fertility	with fertility
Menstruation	✔ with app	✔	✔
Measurement	✔ with app	✔	✔
Temperature reading and cycle day	✔ with app	✔ with app	✔
Fertility indicator		✔	✔
Ovulation display	✔ with app	✔ with app	✔
Volume of stored data	32 days (without syncing)	15+ years	15+ years
Display of past cycle data	with app - all data	with app - all data	360 days
Display of predicted cycle data	80 days	with app - 65 days	6 days
Alarm		✔ with app	✔
Companion app	Teena App	DaysyDay	

Cycle calendar	Teena	Daysy	Lady-Comp
Fertility over past months		✔ with app	✔
Menstruation prediction	with app - approx. 3 cycles	with app - 65 days	on expected. m. days

Cycle and planning statistics	Teena	Daysy	Lady-Comp
Full cycles tracked	✔ with app	✔ with app	✔
Missed measurements		✔ with app	
Average cycle length	✔ with app	✔ with app	✔
Average menstruation length	✔ with app		
Average day of ovulation		✔ with app	✔
Fluctuation range of ovulation in days		✔ with app	✔
Average low temperature		✔ with app	✔
Average high temperature		✔ with app	✔
Average temperature increase after ovulation		✔ with app	✔
Average length of high-temperature phase		✔ with app	✔
Number of monophasic cycles			✔
Number of cycles with CLI (Corpus Luteum Insufficiency)			✔
Planning hits			✔

Pregnancy planning	Teena	Daysy	Lady-Comp
Sexual intercourse input		✔ with app	✔
Possible pregnancy after		15 days	15 days
Probable pregnancy after		18 days	18 days
Date of conception and birth			✔
Pregnancy week			✔

Settings	Teena	Daysy	Lady-Comp
Body Mass Index		✔	✔
Settings °C or °F	✔ with app	✔ with app	✔
Date of birth		✔ with app	✔
Brightness			✔
Volume			✔
Test measurement			✔
Demo mode	✔ for app	✔ for app	✔

Various	Teena	Daysy	Lady-Comp
Powered by	built-in battery	built-in rech. battery	2 x AAA batteries
Life span	5+ years	5+ years	5+ years
Size	4.3 x 1.9 x 0.9 in (H/W/D)	5.1 x 1.9 x 0.9 in (H/W/D)	Ø 3.5 in
Weight	26 g	38 g	130 g